

PHOTOPLETHYSMOGRAPH SIGNAL RECONSTRUCTION BASED ON A NOVEL HYBRID MOTION ARTIFACT DETECTION – REDUCTION APPROACH. PART I: MOTION AND NOISE ARTIFACT DETECTION

PPG and MNA

- Photoplethysmograph (PPG): device to monitor blood volume changes on peripheral tissues. It is used to obtain HR , RR and SpO_2 .
- Motion and Noise Artifacts (MNA): serious obstacles for PPG recording, difficult to filter.
- Sources:
 - *Environmental artifacts*
 - *Physiological artifacts*
 - *Experimental artifacts*
- Causes:
 - *Movement of venous blood and other non-pulsative components along with pulsatile arterial blood*
 - *Variations in the optical coupling between the sensor and the skin*

MNA detection difficulties and how the paper approaches the problem


- They don't have a predetermined frequency band
- Overlapped spectrum with the desire PPG signal
- Use of 4-time domain parameters to quantify MNA as input of a SVM because there are noticeably differences between the values of clean and corrupted parts.
- First step in order to later reconstruct the detected corrupted segments without compromising the clean part of the signal.

Dataset

- Synthetic Dataset for Respiratory Rate Estimation which simulates PPG signals.
- Three respiratory modulations:
 - *Baseline Wander (low frequency artifact in the ECG)*
 - *Amplitude Modulation (proportional to the amplitude of the signal)*
 - *Frequency Modulation (proportional to the frequency of the signal)*
- Sampling Frequency (F_s) = 500 Hz
- Number of samples = 105 000 \rightarrow 210 s
- Data split by respiratory modulation and then divided into training and test phase.

Common procedure for each signal


- Adding gaussian white noise
 - $\text{Noisy Signal} = \text{Original Signal} + \sigma(\text{Signal}) \cdot \text{random value}$
 - To the whole signal with $\text{SNR} = 25 \text{ dB}$
 - To 30 000 samples (60 s), around the 29% of the total, with $\text{SNR} \in [-25 - 10 \text{ dB}]$
- Applying a 3rd order Butterworth bandpass filter with cutoff frequencies of 0.5 and 12 Hz
- Applying a zero phase and reverse filtering
- Splitting the signal into 7-seconds long segments, 30 in total.
- For each segment:
 - Computing 4 parameters
 - Computing labels: 0 = clean, 1 = corrupted (if less than 85% of the segment is free from noise)



Time-Domain Parameters

■ Standard Deviation of Peak-to-Peak Interval


- $STD_{HR} = \sqrt{\frac{1}{N} \sum_{i=1}^N (D_{n,i} - \overline{D_n})^2}$
- $D_{n,i} = T_{peak.n,i} - T_{peak.n,i-1}$
- $\overline{D_n}$ mean peak-to-peak interval of the n th segment



Time-Domain Parameters

■ Standard Deviation of Peak-to-Peak Amplitude


- $STD_{AMP} = \sqrt{\frac{1}{N} \sum_{i=1}^N (A_{n,i} - \overline{A_n})^2}$
- $A_{n,i} = A_{peak.n,i} - A_{trough.n,i+1}$
- $\overline{A_n}$ mean peak-to-peak interval amplitude of the n th segment



Time-Domain Parameters

■ Standard Deviation of Systolic and Diastolic Ratio


- $STD_{SD} = \sqrt{\frac{1}{N} \sum_{i=1}^N (R_{SD,n,i} - \overline{R_{SD,n}})^2}$
- $R_{SD,n,i} = (T_{trough,n,i} - T_{peak,n,i}) / (T_{peak,n,i} - T_{trough,n,i-1})$
- $\overline{R_{SD,n}}$ mean systolic and diastolic time interval ratio of the n th segment




Time-Domain Parameters


■ Standard Deviation of Peak-to-Peak Amplitude

- $STD_{WAV,n} = E[STD_{WAV,n,m}]$
- $STD_{WAV,n,m} = \sqrt{\frac{1}{N} \sum_{i=1}^N (q_{n,i}(m) - \overline{q_{n(m)}})^2}$
- $q_{n,i}(m)$ *mth sample [of N_{samp} taken from a pulse] at the *i*th pulse of the *n*th segment*
- $\overline{q_{n(m)}}$ *mean at the mth sample of the nth segment*




Training Procedure







a,b,c,d) Computed STDs and **e)** Predicted Labels with SNR = -25 dB



a,b,c,d) Computed STDs and **e)** Predicted Labels with SNR = -15 dB



Statistics of STDs values w.r.t. their predicted labels with SNR = -25 dB



Statistics of STDs values w.r.t. their predicted labels with SNR = -15 dB

Parameters obtained from the confusion matrix

Accuracy = (TP + TN)/Tot

Error Rate = (FP + FN)/Tot = 1 – Accuracy

Precision = TP / (TP+FP)


TP = True Positive

TN = True Negative


FP = False Positive

FN = False Negative


Statistics (fm)



Statistics (am)



Statistics (bw)





THE END

THANK YOU FOR LISTENING.